



2015 Release

3B5-2 X3 microprocessor-controlled knee

Reclaim your adventure

The X3 continues to deliver for people looking for adventure, especially if it includes salt (or fresh) water, running, or other prosthetic challenges. As one of the Ottobock family of microprocessor knees, the X3 is the perfect fit for a specific population, eliminating the one size fits all approach to advanced prosthetic fittings.

As a product of a project that pushed the boundaries of prosthetic requirements for active military, the X3 is the answer when “your adventure” includes high mileage or especially rough territory. With new features like the Cockpit app for Android phones, an Activity Report for rehab goal review, and the choice of Intuitive Stance or Deliberate Stance, the X3 delivers more than ever before.

To talk about a demo or trial opportunity, call your Sales Rep at (US) 800 328 4058 or (Canada) 800 665 3327.

New

- ▶ Enhanced Optimized Physiological Gait
- ▶ Choice of two stance functions
- ▶ More control with new Cockpit app
- ▶ Manage rehab with an activity report
- ▶ Quick battery level indicator
- ▶ Free updated X-Soft 1.4

3B5-2 X3 microprocessor-controlled knee

Features and benefits

- A waterproof rating IP68 means X3 is fully submersible, even in salt water
- Up to 5 modes means an easy switch between customized settings to fit unique activities
- Capitalize on improved Optimized Physiological Gait (OPG)
- Confidently descend stairs and step over obstacles
- Choose from Intuitive Stance or Deliberate Stance functions
- Take control with the new Cockpit App for Android devices
- Go farther with up to 5 days of smart battery power management
- Track and report progress with an Activity Report feature
- Make a quick battery check by tipping the knee upside down

Indications

- For bilateral or unilateral patients with knee disarticulation, transfemoral, and hip disarticulation
- For K Level 3 to 4
- Maximum weight is 275 lbs / 125kg
- Contraindicated if the patient's physical and/or mental condition or living situations do not allow the correct handling of X3

Your X3 purchase includes

- X3 knee joint
- 2R19 AXON tube adapter
- Inductive charger
- AC adapter
- Remote control
- Instruction manuals for clinicians (Qualified Personnel) and patients (User)
- Case for charger & AC adapter
- Wallet for cards
- Bluetooth® PIN Card
- Prosthesis Passport
- Cockpit App

Download free X-Soft: <http://ottobock.com/en/datastation>



Technical Data

Suitable for	K3 – K4
Max. body weight	275 lbs / 125 kg

Ordering Information

3B5-2 X3 knee with pyramid adapter
3B5-2=ST X3 knee with screw top

Cover Options

